J. Jimeno

suministros de laboratorio

C/ Vazquez de Menchaca, N°1. Nave 12. Polig. de Argales. 47008 Valladolid.
Tels.: 983 202 342 * 983 202 599. Fax: 983 307 570
www.jjimeno.com - jjimeno@jjimeno.com

FRANCE
CERTIFICAT N°2005101402

Iodolyser

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

TABLE OF CONTENTS

1 – General description of the instrument	p.12
2 – Setting up of the instrument	p.14
3 – Description of the main menus	p.14
4 - Warnings	p.15
5 – Polarisation of the electrode	p.16
6 – Use.....	p.16
6.1. – Free SO2 measure	p.16
6.2. – Total SO2 measure.....	p.17
7 – Notes	p.17
7.1. - Limits of the method of measure.....	p.17
7.2. - Correction to apply in presence of tannic compounds and/or ascorbic acid	p.17
8 – Data transfer.....	p.18
9 – Cleaning and maintenance	p.18
10 – Warranty and return.....	p.18
11 – Technical characteristics	p.19
12 – Declaration of conformity.....	p.19

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A Arcueil, 05/2011

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Instrument :

- | | |
|--|---------------|
| 1. Flask | Ref. 119650/1 |
| 2. GL45 cap with double input | Ref. 119650/2 |
| 3. USB port | |
| 4. Keyboard | |
| 5. Vertical adjustment | |
| 6. Support of electrode and injection | |
| 7. O/I switch | |
| 8. Position of beaker | |
| 9. Screen | |
| 10. Flask/syringe junction silicone hose | Ref. 119650/3 |
| 11. Syringe/injection junction silicone hose | Ref. 119650/3 |
| 12. Position of the platinum electrode | |
| 13. Calibrated injection point | Ref. 119650/4 |
| 14. Electrode tightening screw | Ref. 119650/5 |

Connections :

15. Power cable connector
16. Platinum electrode connector

Accessories :

- | | |
|---|---------------------------|
| 17. Platinum electrode | Ref. 119650/6 |
| 18. Iodide/Iodate KIKIO3 N/64 | Ref. 921500 |
| 19. Sulfuric acid 1/3 | Ref. 941500 |
| 20. Sulfuric acid 1/10 | Ref. 932500 |
| 21. Soda 2N | Ref. 908250 |
| 22. Power cable and adapter | |
| 23. Beakers 50 ml (x3) | Ref. 613050 |
| 24. Graduated pipettes 2 ml, 10 ml, 25 ml | Ref. 604002/604010/604025 |
| 25. Plastic pipette pump | Ref. 606002 |
| 26. Small stirring bars (x2) | Ref. 100062 |

2 – Setting up of the instrument

• Installation of the instrument


- Put the instrument on a flat surface with a plug equipped with a ground connection, complying with current regulations, in a well-ventilated room, isolated from the wine storehouse, protected from direct light and far from a source of heat.
- Fill the glass flask up with the iodide/iodate KIKIO3 N/64
- Put the flask in position in the instrument.

• Installation of the electrode and the injection point

- Remove the protective cap from the electrode and put it on the mobile support.
- Plug the electrode in the right connector
- Use the tightening screw to keep it fixed.

The loops of the electrode must be parallel with the current created by the stirring (cf scheme opposite)

- Insert the injection point into the mobile support.

• Language, date and time setting

- Press **OK** simultaneously with the 0/I switch.
- Choose the language with the drop-down menu and press **OK**.
- Go to the menu: **Time setting**.
- Set up date and time by using $\rightarrow\leftarrow$, and $\downarrow\uparrow$ to increase or decrease numbers.
- Press **OK** and go to the Back menu. Press **OK** to validate.

• Turning on of the instrument

- Plug the instrument on power supply with the adapter and switch it on (0/I switch).
- Place a beaker under the injection point.
- Browse the menus until **Origin**. Press **OK** to initialize the syringe.
- Browse the menus until **Purge** and press **OK**. Repeat the operation once.
- If the injection hose is not completely filled, browse the menus until **Adjustment injection**. Run again the operation **Origin** if required. Run the injection several times in order to fill the hose and the injection point.
- Reinitialize the syringe, using the **Origin** function.

3 – Description of the main menu

• Browsing

- To make the instrument easier to use, drop-down menus are all accessible with the $\uparrow\downarrow$ buttons.
- The button \leftarrow enables to decrease in the menus and so to go back.
- The button **OK** enables to validate or to go into a menu or a submenu.
- Sub-menus also use drop-down menus.

- MENU 1 - **IODO PREPARATION**
- MENU 2 - **ANALYSIS**
 - > Sub-menu 2.1 : **free SO2**
 - > Sub-menu 2.2 : **total SO2**
- MENU 3 - **INJECTION ADJUSTMENT**
- MENU 4 - **AGITATOR**
- MENU 5 - **PURGE**
- MENU 6 - **ORIGINE**
- MENU 7 - **DATA TRANSM.**

- **Detail of the menus**

- MENU 1 – **iodo PREPARATION**

It helps to optimize the performance of the platinum electrode.

- MENU 2 - **ANALYSIS**

It helps to choose the type of analysis to make: free or total SO₂

- > Sub-menu 2.1 : **free SO₂**

- > Sub-menu 2.2 : **total SO₂**

- MENU 3 – **ADJUSTMENT INJECTION**

It enables to adjust the injection before every measurement so that the titrated solution fills the injection hose and the calibrated injection tip.

> TO PROCEED WHEN STARTING THE INSTRUMENT (3)

- MENU 4 - **AGITATOR**

This menu allows to test the good working of motor and to change speed using ↑ and ↓.

- MENU 5 - **PURGE**

This menu enables to purge the syringe and drain the hoses.

> TO PROCEED WHEN STARTING THE INSTRUMENT (2)

- MENU 6 - **ORIGINE**

This menu reinitializes the syringe.

> TO PROCEED WHEN STARTING THE INSTRUMENT (1) AND BEFORE THE FIRST USE.

- MENU 7 – **DATA TRANSM.**

This menu enables to transfer data stocked in the non-volatile memory of the instrument to a computer so as to use such data in an Excel document. The separator of column is the semicolon.

4 - Warnings

- **Use of the instrument**

- Do not place the instrument near a source of heat or directly behind a window.

- In case of problem during a measurement or a polarization, the ↓ button enables to leave the menu and abort the current procedure.

- Purge the syringe when the instrument is left without use for a long time.

- **Solutions**

- The quality and titre of the titration solution KIKIO3 (Ref.921) is very important. Results directly depend on this value.

- Respect the use-by date.

- **Electrode**

- Follow the polarization procedure according to details in chapter 5 before the first use.

- To be efficient, an electrode must have clean parallel round platinum loops.

- The electrode should be rinsed with a wash bottle of distilled water between each titration and dried with clean tissue.

- The platinum loops are malleable. It is possible to easily form them again.

- They can be cleaned with a small brush.

- The electrode must be dry-stored.

5 – Polarisation of the electrode

- **Working condition of the electrode**

To be efficient, an electrode must have clean parallel round platinum loops. It must also be completely polarized.

- **Preparation of the solution for polarisation**

With the supplied pipettes:

- Pour 25 ml of wine or sample into a beaker of 50 ml.
- Add 5 ml of sulfuric acid in 1/3 (Ref. 941).
- Add a stirring bar.

- **Polarization**

- Power the instrument on.
- Put the beaker on the plate of the Iodolyser.
- Adjust the mobile support so that the electrode dips into the solution.
- In the menu **Iodo preparation**, press **OK** to start polarization.
- When the operation is completed, the instrument emits a beep, the polarization is done.

This operation can last several hours.

- **Cleaning**

After polarization, rinse the electrode and the injection point with a wash bottle with distilled water, and dry gently the electrode with a clean tissue.

The platinum loops are malleable. It is possible to easily form them again. They can be cleaned with a small brush.

A polarization is necessary when the instrument has not been used for a long time, and generally once a month.

6 – Use

6.1. Measure of free SO₂

- Prepare a solution of 25 ml of wine or sample to analyze, 5 ml of sulfuric acid in 1/3 (Ref.941) and a stirring bar in a beaker of 50 ml.
- Put the beaker on the plate and adjust the mobile support so that the electrode dips into the solution.
- In the menu **Analysis**, select the sub-menu **Free SO₂** and press **OK**.

The stirring starts and the titration is done automatically. When it is finished, the instrument emits a beep and shows the result on the screen in mg/l.

After titration, rinse the electrode and the injection point with a wash bottle with distilled water, and dry gently the electrode with a clean tissue.

At the end of titration, the instrument reinitializes automatically the syringe.

6.2. Measure of total SO₂

- Prepare a solution of 10 ml of wine or sample to analyze in a beaker of 50 ml.
- Add 2 ml of soda 2N (Ref.908) and wait for 5 minutes.
- In the menu **Analysis**, select the sub-menu **Total SO₂** and press **OK**. The instrument counts down 5 minutes. Press → to pass it.
- Add 20 ml of acid sulfuric in 1/10 (Ref.932) and a stirring bar.
- Put the beaker on the plate and adjust the mobile support so that the electrode dips into the solution.
- Press **OK**.

The stirring starts and the titration is done automatically. When it is finished, the instrument emits a beep and shows the result on the screen in mg/l.

After titration, rinse the electrode and the injection point with a wash bottle with distilled water, and dry gently the electrode with a clean tissue.

At the end of titration, the instrument reinitializes automatically the syringe.

*To be efficient,
an electrode
must have clean
parallel round
platinum loops.*

7 – Notes

7.1. . Limits of the method of measure

This method of measure is used in common practice when the wines to analyze are enough far from their legal limit in SO₂.

Otherwise, and for a better accuracy, it is advised to use the official method, as recommended by the OIV, which is the distillation-oxidation with the Frantz Paul equipment.

7.2. Correction to apply in presence of tannic compounds and ascorbic acid

In case of highly colored or very tannic red wines, and/or wines with ascorbic acid, it is suitable to make a correction on the Free SO₂, by doing the following operation:

- Prepare a solution of 25 ml of wine or sample to analyze in a beaker of 50 ml. Add 5 ml of solution of propionaldehyde (Ref.973) (non supplied) and wait for 15 minutes.
- Add 5 ml of acid sulfuric in 1/3 (Ref.941) and a stirring bar.
- Put the beaker on the plate and adjust the mobile support so that the electrode dips into the solution.
- Use the titration mode **ANALYSIS > Free SO₂**

The stirring starts and the titration is done automatically. When it is finished, the instrument emits a beep and shows the result on the screen in mg/l.

Deduce this value from the one you obtain when doing the simple free SO₂ analysis to obtain the corrected concentration of free SO₂.

After titration, rinse the electrode and the injection point with a wash bottle with distilled water, and dry gently the electrode with a clean tissue.

At the end of titration, the instrument reinitializes automatically the syringe.

8 – Data transfer

This menu enables to transfer data stocked in the non-volatile memory of the instrument to a computer with a USB port.

- Connect a series cable between the instrument and your computer
- Transfer data with the programme.
- Open Excel and open the transferred file. Precise that the separator of column is the semi colon.

La USB connection can also be used to update the programme used by the instrument.

9 – Cleaning and maintenance

- **Cleaning of the injection system**

When the instrument is not used for more than 8 days : remove the flask, empty it or save the solution. Place a beaker with distilled water and purge twice with the menu Purge.

When the instrument is not used for more than 1 month : remove the flask, empty it or save the solution. Fill the flask with distilled water and purge 5 times. Then remove the flask and purge twice to drain the hoses.

- **Cleaning of the electrode**

See chapter 4 «Warnings».

10 – Warranty and return

- **Instruments under warranty should be returned:**

- In their original packaging (parcels and styrofoam);
- With shipping costs paid, including insurance against any breakage during the carriage.

We do not warranty any failure due to an ordinary wear, a non-approved use, an on-site service on the instruments done by a third After Sales Service or to an external accident.

- **Instruments out of warranty should be returned:**

- In their original packaging when possible or perfectly secured in a double packing to avoid any direct contact between the different elements;
- With shipping costs paid, including insurance against any breakage during the carriage.

- **The warranty applies to :**

- All parts used normally, except parts of everyday use such as glassworks, electrode, leds, switches and connectors, along with consumables.

11 – Technical characteristics

Iodolyser - Ref.119650

Type of sample	Wines, lees for distillations, fruit juices, alcoholic drinks, vinegars
Programmes	Free SO ₂ , Total SO ₂ , Correction to apply in presence of tannic compounds and ascorbic acid, Polarisation of the electrode
Measure scale	mg/L of SO ₂
Resolution	1 mg/L
Capacity of syringe	10 ml
Precision of seringue	10 µl per injection
Electrode	Double platinum electrode
Valve	Double pinch valve
Keyboard	5 keys
Screen	LCD with backlighting
Languages	French, English, Spanish, Italian, German
Input/output	USB (to download data on PC and upload upgrade of the software on Iodolyser)
Dimensions and weight	24 x 35 x 40 cm - 4 kg
Power supply	110-220V - 50/60 Htz
Adapter	15V

12 – Declaration of conformity

LABORATOIRES DUJARDIN-SALLERON

18, rue Henri Barbusse

94117 ARCUEIL CEDEX

FRANCE

Tel : 00 33 1-45-46-04-05

Fax : 00 33 1-45-46-01-13

E-mail : info@dujardin-salleron.com

STATEMENT OF CONFORMITY 

Type of instrument: Apparatus for SO₂ determination

Model IODOLYSER Ref.119650

Has been developed, designed and produced according to the regulations of the Directive CEM 2004/108/CE and of the Directive Basse tension 2006/95/CE.

In accordance with the harmonized standards:

EN 55022:2006	IEC 61000-4-4:2004
EN 61000-3-2:2006	IEC 61000-4-11:2004
EN 61000-3-3:1995+A1:2001+A2:2005	IEC 61000-4-3:2006
EN 55024:1998+A1:2001+A2:2003	IEC 61000-4-5:2005
EN 61010-1:2004	

Arcueil, 05/2011
Laurent Dubreuil, President